

Tissue isolation and embedding in OCT (Optimal Cutting Temperature Compound)

Lessard Group, GUDMAP Consortium

1. Dissect tissue (e.g. bladder, ureter, etc.) rapidly and place each specimen in OCT using a plastic mold (Fisher Scientific).
2. Be sure the specimen is covered in central position, not too near the top and in the desired orientation for sectioning.
3. Immediately freeze the sample by placing it on dry ice.
4. When the OCT is completely frozen keep the mold in dry ice, and then store long term in liquid nitrogen freezer or at -80°C .