

Processing slides for LCM

Potter Group, GUDMAP Consortium

1. Remove slides from -80°C freezer, store temporarily on dry ice, air dry on metal slide warmer turned off or set to about 30 to 35°C , for 3 min.
2. Fix two minutes. 1:1 mix of acetone 75% ethanol gives a nice compromise, with good section sticking and reasonable histology, and good lectin staining.
3. Rinse 2 min in 1/10X PBS, ice cold, to dissolve off some OCT. The 1/10 PBS is made by diluting 1X PBS ten fold with sterile autoclaved super water.
4. Stain about 6 min in lectin, on ice. For PNA use 5 μl per ml of 1/10 PBS. Stain on ice.
5. Rinse 1/10 PBS, ice cold. Two gentle 10 sec dips and then 3 min.
6. Dehydrate, in ethanol, 75%, 95%, 100%, 10-15 sec each, with 2-3 dips, and then another 100% one minute.
7. Xylene 1.5 min with three dips.
8. Xylene 2 min with two dips.
9. Air dry about 2-4 min.