

E18.5 mouse DRG dissociation for FACS, 060214

Masato Hoshi, MD, PhD,
Sanjay Jain Lab, Renal Division, Department of Medicine,
Washington University School of Medicine

Reagents;

- Hibernate-E; Life Technologies, #A12476-01
- DMEM/F12; Life Technologies, #11330
- 0.05% Trypsin/0.02% EDTA; Sigma, #59417C or #T3924
- DNaseI; Sigma, #D5025
- Collagenase; Sigma, #C0130
- RNase inhibitor; RNaseOUT (Life Technologies, 40units/ul) or equal
- 7-AAD; Life Technologies, #A1310 (powder, dissolve in DMSO at 1mg/ml)

Procedures;

- 1) Euthanize a pregnant mouse following your animal protocol approval and extract embryos from the uterus.
- 2) Save the embryos/fetuses in a 12-well plate with 3ml of Hibernate-E medium (Life Technologies) in each well of the plate on ice.
- 3) Take one embryo in 10cm plastic plate with dissection medium (DMEM/F12 (Life Technologies) or Hibernate-E (30ml)) and extirpate DRGs from embryos. (Th11 to S1, total 20 DRGs from each embryo; or as needed for your study)
(*when you use DMEM/F12, the medium should be put in a CO₂ incubator (cell culture, 5% CO₂ incubator at 37C is fine) for more than 2 hours to keep the pH balanced (gauge by the color, prefer orange).
- 4) Collect DRGs in a 1.5ml tube on ice with dissection medium in a pipette tip. Keep the tube on ice until all dissections are done.
- 5) Centrifuge at 500 xg for 5min at 4C and discard as much supernatant as possible by using a pipette.
- 6) Add 500ul (to 1ml, depending on the number of DRGs) of 0.05% Trypsin/ 0.02% EDTA (Sigma) with 200ug/ml of DNaseI (Sigma) to the tube and briefly vortex to mix.
- 7) Incubate it at 37C for 15min.
- 8) Triturate with P-1000, then P-200 pipette and dissociate the specimens. Still see some chunks.
- 9) Add 500ul of DMEM/F12 + 10% FCS (filtered) and mix well by pipetting.
- 10) Centrifuge at 500 xg for 5min at 4C and discard supernatant by using a pipette.
- 11) Add 400ul of DMEM/F12 + 10% FCS (filtered), mix well by pipetting, and then add 600ul of collagenase solution (Sigma, dissolved in DMEM/F12 + 10%FCS (sterilized) at 2mg/ml; final collagenase concentration is 1.2mg/ml). Mix well by pipetting.
- 12) Incubate it at 37C for 10min.
- 13) Triturate with P-200 pipette till you notice complete dissociation of the specimen.
- 14) Centrifuge at 500 xg for 5min at 4C and discard supernatant by using a pipette. Leave around 30ul to avoid disturbing the cell pellet.

- 15) Add 500ul of PBS + 5% FCS (filtered) with 0.1units/ul of RNase inhibitor (Life Technologies, 40units/ul, 25ul to 10ml PBS/FCS solution) and mix well by pipetting.
- 16) Centrifuge at 500 xg for 5min at 4C and discard supernatant by using a pipette.
- 17) Repeat this washing step twice more with PBS/FCS/RNase inhibitor solution. Total three times. In the last wash, the cells get dispersed very quickly with only few rounds of pipetting.
- 18) After discarding as much supernatant as possible, add 300ul (or appropriate volume for FACS) of PBS + 5% FCS (filtered) with 0.1units/ul of RNase inhibitor and mix well by pipetting.
- 19) Filter the dissociated cells with 40um pore size cell strainer (BD, #352340) and collect the cells into a FACS tube (BD, #352063).
- 20) Add 7-AAD at 1ug/ml concentration, if needed, and do sorting. 7-AAD is added 10min before sorting. We have had good success using MoFlo, Sony or Ariall, at 100 um nozzle size, and 30psi.
- 21) Sorted cells are collected in 750ul of TRIzol-LS (Life Technologies, #10296-010) for RNA extraction and saved in -80C until use.