

# Library construction protocol for nanogram amounts of total RNA for RNAseq – HIGH VOLUME

## Jain Group Protocols, GUDMAP Consortium

1) 3ng, 10ng or 30ng total RNA (in 9µl) was reverse transcribed and amplified using the SMARTer® Ultra™ Low Input RNA for Illumina® Sequencing - HV (Clontech Laboratories, cat# 634828) according to the manufacturer's instructions.

2) The cDNA was sheared for 3 minutes with the Covaris ultra-sonicator.

3) The libraries were prepared using VeraSeq™ ULtra DNA Polymerase (Enzymatics, cat# P7520L) for 12 cycles. PCR conditions were 98°C for 30s; 12 cycles of 98°C for 10s, 64°C for 30s, 72°C for 30s; 72°C for 5 min and then hold at 4°C. The forward primer used was: AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT

The reverse primer with the indexes used were:

| Reverse Primer                                                                 | Reverse Complement | Reported index & Sample Name |
|--------------------------------------------------------------------------------|--------------------|------------------------------|
| 64,<br>CAAGCAGAAGACGGCATACGAGAT <b>GTAATCA</b> GTGACTGGAGTTCAGACGTGTGCTCTTCCGA | TGATTAC            | TGATTAC-30ng_High-Volume     |
| 65,<br>CAAGCAGAAGACGGCATACGAGAT <b>CCCAGCA</b> GTGACTGGAGTTCAGACGTGTGCTCTTCCGA | TGCTGGG            | TGCTGGG-10ng_High-Volume     |
| 66,<br>CAAGCAGAAGACGGCATACGAGAT <b>CTGTGTC</b> GTGACTGGAGTTCAGACGTGTGCTCTTCCGA | GACACAG            | GACACAG-3ng_High-Volume      |

4) The purified DNA yield was 16.5ng/µl in 30µl and sequencing was done on the Illumina HiSeq2500 following the manufacturer's protocols.