

Fluorescent Immunocytochemistry Protocol on Frozen sections

Gaido Group, GUDMAP Consortium

1. Prepare cryostat section at 5-7 mm (for In situ, may 10 mm)
2. Fix with 4% paraformaldehyde (diluted 16% paraformaldehyde at 1:4 in PBS), 10 min, room temperature
3. Wash with PBS/Triton, 5 min X 3
4. Block with Blocking buffer (500 ml normal goat serum to 10 ml PBS + 10 ml Triton, serum match to species the second antibody from), 1 hr, room temperature
5. Wash with PBS, 5 min X 1
6. Primary antibody diluted in 3% BSA (300mg in 10 ml PBS), overnight, at 4 °C
7. Wash with PBS, 5 min X 3
8. Incubate with biotinylated secondary antibody, 1:500 dilution (pay attention to different species, against primary), 1 hr, room temperature
9. Wash with PBS, 5 min X 3
10. Mount coverslips with VectaMount
11. Visualization using microscope
12. Analysis using Confocal Software